


PROBE

Ø 2.0mm

20.0cm (8")



SKU: 080-004-20

Categories: [Open Surgery](#), [Probes](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Buttoned Probes](#)

Aesculap BN140

KLS Martin 16-100-20

Medicon 09.05.20

SPECIFICATIONS [without eye](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Buttoned Probes](#)

Aesculap BN140

KLS Martin 16-100-20

Medicon 09.05.20

SPECIFICATIONS

[without eye](#)