


PROBE

Ø 1.0mm

16.0cm (6.4")



SKU: 080-000-16

Categories: [Open Surgery](#), [Probes](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Buttoned Probes
Aesculap	BN116
KLS Martin	16-102-16
Medicon	09.04.16
SPECIFICATIONS	without eye

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Buttoned Probes
Aesculap	BN116
KLS Martin	16-102-16
Medicon	09.04.16

SPECIFICATIONS

[without eye](#)