

PROBE WITH EYE

Ø 2.0mm

16.0cm (6.4")

SKU: 080-013-16

Categories: [Open Surgery](#), [Probes](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Buttoned Probes
Aesculap	BN036
KLS Martin	16-108-16
Medicon	09.10.16
SPECIFICATIONS	with eye

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Buttoned Probes
Aesculap	BN036
KLS Martin	16-108-16
Medicon	09.10.16

SPECIFICATIONS

[with_eye](#)