

## PROBE WITH EYE

Ø 2.0mm

20.0cm (8")

**SKU:** 080-013-20

**Categories:** [Open Surgery](#), [Probes](#)



## ZUSÄTZLICHE INFORMATIONEN

|                       |                                 |
|-----------------------|---------------------------------|
| <b>UNIT</b>           | <a href="#">pcs.</a>            |
| <b>TYPE</b>           | <a href="#">Buttomed Probes</a> |
| <b>Aesculap</b>       | BN040                           |
| <b>KLS Martin</b>     | 16-108-20                       |
| <b>Medicon</b>        | 09.10.20                        |
| <b>SPECIFICATIONS</b> | <a href="#">with eye</a>        |

## ZUSÄTZLICHE INFORMATIONEN

|                   |                                 |
|-------------------|---------------------------------|
| <b>UNIT</b>       | <a href="#">pcs.</a>            |
| <b>TYPE</b>       | <a href="#">Buttomed Probes</a> |
| <b>Aesculap</b>   | BN040                           |
| <b>KLS Martin</b> | 16-108-20                       |
| <b>Medicon</b>    | 09.10.20                        |

## SPECIFICATIONS

[with\\_eye](#)