

PROBE WITH EYE

Ø 2.0mm

20.0cm (8")

SKU: 080-013-20

Categories: [Open Surgery](#), [Probes](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Buttoned Probes
Aesculap	BN040
KLS Martin	16-108-20
Medicon	09.10.20
SPECIFICATIONS	with eye

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Buttoned Probes
Aesculap	BN040
KLS Martin	16-108-20
Medicon	09.10.20

SPECIFICATIONS

[with_eye](#)