

PROBE WITH EYE

Ø 2.0mm

11.5cm (4.6")

SKU: 080-013-11

Categories: [Open Surgery](#), [Probes](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Buttoned Probes](#)

KLS Martin 16-108-11

SPECIFICATIONS [with eye](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Buttoned Probes](#)

KLS Martin 16-108-11

SPECIFICATIONS [with eye](#)