

PUNCH; THROUGH-CUTTING; ALLIGATOR; Ø 1.6MM

Fr. 5; 34cm; semi-rigid; single action

Luer-Lock

SKU: 583-258-16LL

Categories: [Endoscopy](#), [GUBBINI/Gyn](#), [Gynecology](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
RANGE [Micro Instruments Hysteroscopy](#)
TYPE ["Semi-Rigid" with LuerLock](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
RANGE [Micro Instruments Hysteroscopy](#)
TYPE ["Semi-Rigid" with LuerLock](#)