


## PUSHER FOR ENDOPROTHESIS

with goldring, x-ray

2.3mm, length 270cm; single use



**SKU:** PK-FPU7270

**Category:** [Endoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Instruments for Double-Balloon-Enteroscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Instruments for Double-Balloon-Enteroscopy](#)