

RAMPLEY SPONGE FORCEPS

straight; serrated

18.0cm (7.2")


SKU: 030-190-18

Categories: [Hemostats](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---|
| UNIT | pcs. |
| TYPE | Sponge and Dressing Forceps |
| Aesculap | BF116R |
| AUTHORS | RAMPLEY |
| SPECIFICATIONS | serrated , straight |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---|
| UNIT | pcs. |
| TYPE | Sponge and Dressing Forceps |
| Aesculap | BF116R |
| AUTHORS | RAMPLEY |
| SPECIFICATIONS | serrated , straight |