

## RESECTOSCOPY SHEATH

13Fr.; 1 turnable stopcock

pediatric



**SKU:** 630-510-13

**Categories:** [Endoscopy](#), [Urology](#)

### ZUSÄTZLICHE INFORMATIONEN

|              |  |
|--------------|--|
| <b>UNIT</b>  | <a href="#">pcs.</a>                   |
| <b>RANGE</b> | <a href="#">Pediatric Resectoscopy</a> |
| <b>TYPE</b>  | <a href="#">Sheaths</a>                |

### ZUSÄTZLICHE INFORMATIONEN

|              |  |
|--------------|--|
| <b>UNIT</b>  | <a href="#">pcs.</a>                   |
| <b>RANGE</b> | <a href="#">Pediatric Resectoscopy</a> |
| <b>TYPE</b>  | <a href="#">Sheaths</a>                |