

## RETRIEVAL BASKET, STAINLESS

hexagonal; 20 X 60mm; 4 wires

2.5mm.; 210cm; Lithotripsy



**SKU:** PK-1SHPK2520210AK

**Category:** [Endoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Stonebaskets for Retrieval & Lithotripsy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Stonebaskets for Retrieval & Lithotripsy](#)