

## RETRIEVAL BASKET, STAINLESS

hexagonal; 30 x 60; 6 wires

2.5mm.; 210cm; Lithotripsy



**SKU:** PK-1SHPK2530210VAK

**Category:** [Endoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Stonebaskets for Retrieval & Lithotripsy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Stonebaskets for Retrieval & Lithotripsy](#)