


RETRIEVAL BASKET, STAINLESS

spiral; 20 X 40mm; 4 wires

2.5mm; 210cm, Lithotripsy



SKU: PK-1SSPK2520210AK

Category: [Endoscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Stonebaskets for Retrieval & Lithotripsy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Stonebaskets for Retrieval & Lithotripsy](#)