

## RETRIEVAL BASKET, STAINLESS

hexagonal, 30 x 60mm, 4 wires

2.3mm; 230cm



**SKU:** PK-1BHPK2330230K

**Category:** [Endoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Stonebaskets for Retrieval & Lithotripsy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Stonebaskets for Retrieval & Lithotripsy](#)