

RETRIEVAL BASKET, STAINLESS

spiral, 30 x 60mm, 4 wires

2.3mm; 230cm



SKU: PK-1BSPK2330230K

Category: [Endoscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Stonebaskets for Retrieval & Lithotripsy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Stonebaskets for Retrieval & Lithotripsy](#)