

RETRIEVAL BASKET, STAINLESS STEEL

hexagonal, 20 x 40mm, 4 wires, hexagonal

2.3mm.; 230cm (single use)



SKU: PK-1BHEW2320230K-S

Category: [Endoscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Single Use Instruments](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Single Use Instruments](#)