

## RETRIEVAL BASKET, STAINLESS STEEL

hexagonal, 20 x 40mm, 4 wires, hexagonal

2.3mm.; 230cm (single use)



**SKU:** PK-1BHEW2320230K-S

**Category:** [Endoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Single Use Instruments](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Single Use Instruments](#)