

## RING-APPLICATOR FOR OP LAPAROSCOPE

Ø 6.0mm; 335mm

**SKU:** 587-235-30

**Categories:** [Endoscopy](#), [Gynecology](#)



### ZUSÄTZLICHE INFORMATIONEN

|              |   |
|--------------|---|
| <b>UNIT</b>  | <a href="#">pcs.</a>                    |
| <b>RANGE</b> | <a href="#">Laparoscopy Instruments</a> |
| <b>TYPE</b>  | <a href="#">Ring Aplicators</a>         |

### ZUSÄTZLICHE INFORMATIONEN

|              |   |
|--------------|---|
| <b>UNIT</b>  | <a href="#">pcs.</a>                    |
| <b>RANGE</b> | <a href="#">Laparoscopy Instruments</a> |
| <b>TYPE</b>  | <a href="#">Ring Aplicators</a>         |