

## RING-APPLICATOR FOR OP LAPAROSCOPE

Ø 6.0mm; 335mm

**SKU:** 587-235-30

**Categories:** [Endoscopy](#), [Gynecology](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Laparoscopy Instruments</a>
<b>TYPE</b>	<a href="#">Ring Aplicators</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Laparoscopy Instruments</a>
<b>TYPE</b>	<a href="#">Ring Aplicators</a>