

## RING STRIPPER, ROUND MODEL

Ø 3.0mm

50.0cm (20")

**SKU:** 080-420-03

**Categories:** [Open Surgery](#), [Probes](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Vascular Strippers</a>
<b>KLS Martin</b>	16-251-03

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Vascular Strippers</a>
<b>KLS Martin</b>	16-251-03