

## ROCHESTER-PEAN HAEMOSTATIC FORCEPS

straight; serrated

24.0cm (9.6")



**SKU:** 040-260-24

**Categories:** [Hemostats](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Artery Forceps</a>
<b>Aesculap</b>	BH450
<b>KLS Martin</b>	13-410-24
<b>AUTHORS</b>	<a href="#">PEAN</a> , <a href="#">ROCHESTER</a>
<b>Medicon</b>	16.10.24
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Artery Forceps</a>
<b>Aesculap</b>	BH450
<b>KLS Martin</b>	13-410-24

**AUTHORS**

[PEAN, ROCHESTER](#)

**Medicon**

16.10.24

**SPECIFICATIONS**

[serrated](#), [straight](#)