

ROCHESTER-PEAN HAEMOSTATIC FORCEPS

straight; serrated

24.0cm (9.6")

SKU: 040-260-24

Categories: [Hemostats](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps
Aesculap	BH450
KLS Martin	13-410-24
AUTHORS	PEAN , ROCHESTER
Medicon	16.10.24
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps
Aesculap	BH450
KLS Martin	13-410-24

AUTHORS

[PEAN, ROCHESTER](#)

Medicon

16.10.24

SPECIFICATIONS

[serrated](#), [straight](#)