

## ROCHESTER-PEAN HAEMOSTATIC FORCEPS

straight; serrated

30.0cm (12")



**SKU:** 040-260-30

**Categories:** [Hemostats](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Artery Forceps</a>
<b>Aesculap</b>	BH474
<b>KLS Martin</b>	13-410-30
<b>AUTHORS</b>	<a href="#">PEAN</a> , <a href="#">ROCHESTER</a>
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Artery Forceps</a>
<b>Aesculap</b>	BH474
<b>KLS Martin</b>	13-410-30
<b>AUTHORS</b>	<a href="#">PEAN</a> , <a href="#">ROCHESTER</a>

## SPECIFICATIONS

[serrated](#), [straight](#)