

ROGERS HYSTERECTOMY FORCEPS

straight; serrated

25.5cm (10")



SKU: 040-580-24

Categories: [Hemostats](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Hysterectomy Forceps
Aesculap	BJ550
KLS Martin	13-637-10
AUTHORS	ROGERS
Medicon	16.35.50
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Hysterectomy Forceps
Aesculap	BJ550
KLS Martin	13-637-10

AUTHORS

[ROGERS](#)

Medicon

16.35.50

SPECIFICATIONS

[serrated](#), [straight](#)