

ROGERS HYSTERECTOMY FORCEPS

angled; serrated

21.0cm (8 1/2")




SKU: 040-585-21

Categories: [Hemostats](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Hysterectomy Forceps
Aesculap	BJ533
KLS Martin	13-635-40
AUTHORS	ROGERS
SPECIFICATIONS	angled , serrated



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Hysterectomy Forceps
Aesculap	BJ533
KLS Martin	13-635-40
AUTHORS	ROGERS

SPECIFICATIONS

[angled](#), [serrated](#)