

## ROGERS HYSTERECTOMY FORCEPS

angled; serrated

21.0cm (8 1/2")




**SKU:** 040-585-21

**Categories:** [Hemostats](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Hysterectomy Forceps</a>
<b>Aesculap</b>	BJ533
<b>KLS Martin</b>	13-635-40
<b>AUTHORS</b>	<a href="#">ROGERS</a>
<b>SPECIFICATIONS</b>	<a href="#">angled</a> , <a href="#">serrated</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Hysterectomy Forceps</a>
<b>Aesculap</b>	BJ533
<b>KLS Martin</b>	13-635-40
<b>AUTHORS</b>	<a href="#">ROGERS</a>

## SPECIFICATIONS

[angled](#), [serrated](#)