

ROGERS HYSTERECTOMY FORCEPS

angled; serrated

24.0cm (9.6")


SKU: 040-585-24

Categories: [Hemostats](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Hysterectomy Forceps
Aesculap	BJ553
KLS Martin	13-637-40
AUTHORS	ROGERS
Medicon	16.35.53
SPECIFICATIONS	angled , serrated



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Hysterectomy Forceps
Aesculap	BJ553
KLS Martin	13-637-40

AUTHORS

[ROGERS](#)

Medicon

16.35.53

SPECIFICATIONS

[angled](#), [serrated](#)