

## ROGERS HYSTERECTOMY FORCEPS

angled; serrated

24.0cm (9.6")




**SKU:** 040-585-24

**Categories:** [Hemostats](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Hysterectomy Forceps</a>
<b>Aesculap</b>	BJ553
<b>KLS Martin</b>	13-637-40
<b>AUTHORS</b>	<a href="#">ROGERS</a>
<b>Medicon</b>	16.35.53
<b>SPECIFICATIONS</b>	<a href="#">angled</a> , <a href="#">serrated</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Hysterectomy Forceps</a>
<b>Aesculap</b>	BJ553
<b>KLS Martin</b>	13-637-40

## AUTHORS

[ROGERS](#)

## Medicon

16.35.53

## SPECIFICATIONS

[angled](#), [serrated](#)