

ROGERS HYSTERECTOMY FORCEPS

straight; serrated

21.5cm (8.6")




SKU: 040-580-21

Categories: [Hemostats](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Hysterectomy Forceps
Aesculap	BJ530
KLS Martin	13-635-10
AUTHORS	ROGERS
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Hysterectomy Forceps
Aesculap	BJ530
KLS Martin	13-635-10
AUTHORS	ROGERS

SPECIFICATIONS

[serrated](#), [straight](#)