

## ROSEN ROUND CUTTING KNIFE

Ø 3.0mm; 45° angled;

16.0cm (6.4")


**SKU:** 111-666-30

**Categories:** [Open Surgery](#), [Otology](#)



### ZUSÄTZLICHE INFORMATIONEN

|                   |                        |
|-------------------|------------------------|
| <b>UNIT</b>       | <a href="#">pcs.</a>   |
| <b>TYPE</b>       | <a href="#">Knives</a> |
| <b>Aesculap</b>   | OG132                  |
| <b>KLS Martin</b> | 36-022-01              |
| <b>AUTHORS</b>    | <a href="#">ROSEN</a>  |



### ZUSÄTZLICHE INFORMATIONEN

|                   |                        |
|-------------------|------------------------|
| <b>UNIT</b>       | <a href="#">pcs.</a>   |
| <b>TYPE</b>       | <a href="#">Knives</a> |
| <b>Aesculap</b>   | OG132                  |
| <b>KLS Martin</b> | 36-022-01              |
| <b>AUTHORS</b>    | <a href="#">ROSEN</a>  |