

SANTULLI BULLDOG CLAMP

angled; atraumatic

7.0cm (2.8")


SKU: 026-195-07

Categories: [Hemostats](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bulldog Clamps
KLS Martin	24-295-07
AUTHORS	SANTULLI
SPECIFICATIONS	angled , atraumatic



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bulldog Clamps
KLS Martin	24-295-07
AUTHORS	SANTULLI
SPECIFICATIONS	angled , atraumatic