

## SAROT HAEMOSTATIC FORCEPS

straight; serrated

24.0cm (9.6")



**SKU:** 040-450-24

**Categories:** [Hemostats](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Artery Forceps</a> , <a href="#">Bronchus Forceps</a>
<b>Aesculap</b>	BH224
<b>KLS Martin</b>	13-432-24
<b>AUTHORS</b>	<a href="#">SAROT</a>
<b>Medicon</b>	16.14.74
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Artery Forceps</a> , <a href="#">Bronchus Forceps</a>
<b>Aesculap</b>	BH224
<b>KLS Martin</b>	13-432-24

## AUTHORS

[SAROT](#)

## Medicon

16.14.74

## SPECIFICATIONS

[serrated](#), [straight](#)