

## SCHMID IRRIGATION CANNULA

Ø 3.0mm; rigid

15.0cm (6")

**SKU:** 086-140-03

**Categories:** [Open Surgery](#), [Puncture / Special Needles](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Special Needles</a>
<b>Aesculap</b>	FB153
<b>KLS Martin</b>	18-695-03
<b>AUTHORS</b>	<a href="#">SCHMID</a>
<b>Medicon</b>	09.96.53

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Special Needles</a>
<b>Aesculap</b>	FB153
<b>KLS Martin</b>	18-695-03
<b>AUTHORS</b>	<a href="#">SCHMID</a>

**Medicon**

09.96.53