

SCHMIEDEN-PAYR LIGATURE CONDUCTOR

curved

19.0cm (7.6")


SKU: 056-506-19

Categories: [Open Surgery](#), [Suture](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Ligature Conductors
Aesculap	BM905
KLS Martin	20-871-19
AUTHORS	PAYR , SCHMIEDEN
Medicon	13.32.19
SPECIFICATIONS	curved



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Ligature Conductors
Aesculap	BM905
KLS Martin	20-871-19

AUTHORS

[PAYR](#), [SCHMIEDEN](#)

Medicon

13.32.19

SPECIFICATIONS

[curved](#)