

SCHNIDT TONSIL FORCEPS

curved; open ring

19.0cm (7.6")



SKU: 132-187-19

Categories: [Open Surgery](#), [Tonsils](#), [Tracheotomy](#), [Tongue](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tonsils Forceps
Aesculap	BH961
KLS Martin	39-447-19
AUTHORS	SCHNIDT
Medicon	17.21.50



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tonsils Forceps
Aesculap	BH961
KLS Martin	39-447-19
AUTHORS	SCHNIDT

Medicon

17.21.50