

SCHUKNECHT SUCTION TUBE

20 Gauge (0.9mm); WL: 75mm

SKU: 086-310-09

Categories: [Open Surgery](#), [Puncture / Special Needles](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Suction Tubes
AUTHORS	SCHUKNECHT

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Suction Tubes
AUTHORS	SCHUKNECHT