

SCHUKNECHT SUCTION TUBE

20 Gauge (0.9mm); WL: 75mm



SKU: 086-310-09

Categories: [Open Surgery](#), [Puncture / Special Needles](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Suction Tubes
AUTHORS	SCHUKNECHT

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Suction Tubes
AUTHORS	SCHUKNECHT