

SCISSORS; POINTED; Ø 1.6MM; FR. 5; 34CM

semi-rigid; single action; Luer-Lock

with handle

SKU: 583-229-16LL

Categories: [Endoscopy](#), [GUBBINI/Gyn](#), [Gynecology](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Micro Instruments Hysteroscopy
TYPE	"Semi-Rigid" with LuerLock



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Micro Instruments Hysteroscopy
TYPE	"Semi-Rigid" with LuerLock