

## SEALING CAP RED

for trocar sleeve Ø 5.5mm

10 pcs./box; Price per box



**SKU:** 626-170-50

**Categories:** [Endoscopy](#), [Laparoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">PU</a>
<b>RANGE</b>	<a href="#">Trocar Systems</a>
<b>TYPE</b>	<a href="#">Sheaths</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">PU</a>
<b>RANGE</b>	<a href="#">Trocar Systems</a>
<b>TYPE</b>	<a href="#">Sheaths</a>