

SEALING SET Ø 20.0MM; CONSISTING OF:

sealing cap grey; inner sealing grey

O-ring grey; 10 pcs./box; Price per box



SKU: 626-171-20

Categories: [Endoscopy](#), [Laparoscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	PU
RANGE	Trocar Systems
TYPE	Sheaths

ZUSÄTZLICHE INFORMATIONEN

UNIT	PU
RANGE	Trocar Systems
TYPE	Sheaths