

SEIFFERT EXCHANGEABLE TIP

toothed

for BRÜNINGS interchangeable tubes


SKU: 131-482-00

Categories: [Goitre](#), [Laryngology](#), [Mouth](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Universal handles and tips
Aesculap	OP202
KLS Martin	39-604-01
AUTHORS	SEIFFERT



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Universal handles and tips
Aesculap	OP202
KLS Martin	39-604-01
AUTHORS	SEIFFERT