

SEMKEN DRESSING FORCEPS

straight; serrated

12.5cm (5"); width of jaw 1.3 +/-0.2

SKU: 020-040-12

Categories: [Forceps](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---|
| UNIT | pcs. |
| TYPE | Dressing Forceps , Tissue Forceps |
| Aesculap | BD215 |
| KLS Martin | 12-172-12 |
| AUTHORS | SEMKEN |
| Medicon | 06.30.12 |
| SPECIFICATIONS | serrated , straight |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|---|
| UNIT | pcs. |
| TYPE | Dressing Forceps , Tissue Forceps |
| Aesculap | BD215 |
| KLS Martin | 12-172-12 |

AUTHORS

[SEMKEN](#)

Medicon

06.30.12

SPECIFICATIONS

[serrated](#), [straight](#)