


## SILICON RING

autoclavable; medium

27x3mm



**SKU:** 196-988-27

**Categories:** [Open Surgery](#), [Orthopaedics](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Lead Hands</a>
<b>KLS Martin</b>	23-960-27
<b>Medicon</b>	78.000.06

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Lead Hands</a>
<b>KLS Martin</b>	23-960-27
<b>Medicon</b>	78.000.06