


## SILICON RING

autoclavable; small

24x3mm



**SKU:** 196-988-24

**Categories:** [Open Surgery](#), [Orthopaedics](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Lead Hands</a>
<b>KLS Martin</b>	23-960-24
<b>Medicon</b>	78.00.05

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Lead Hands</a>
<b>KLS Martin</b>	23-960-24
<b>Medicon</b>	78.00.05