

SNARE WIRE FOR TONSIL SNARE

Ø 0.3mm; package of 12 pieces

SKU: 132-521-03

Categories: [Open Surgery](#), [Tonsils](#), [Tracheotomy](#), [Tongue](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tonsil Snares
Aesculap	OM750
KLS Martin	39-513-01

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tonsil Snares
Aesculap	OM750
KLS Martin	39-513-01