

## SNARE WIRE FOR TONSIL SNARE

Ø 0.3mm; package of 12 pieces

**SKU:** 132-521-03

**Categories:** [Open Surgery](#), [Tonsils](#), [Tracheotomy](#), [Tongue](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Tonsil Snares</a>
<b>Aesculap</b>	OM750
<b>KLS Martin</b>	39-513-01

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Tonsil Snares</a>
<b>Aesculap</b>	OM750
<b>KLS Martin</b>	39-513-01