

## SPENCER-WELLS HAEMOSTATIC FORCEPS

straight; serrated

14.0cm (5.6")



**SKU:** 040-290-14

**Categories:** [Hemostats](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Artery Forceps</a>
<b>KLS Martin</b>	13-248-14
<b>AUTHORS</b>	<a href="#">SPENCER</a> , <a href="#">WELLS</a>
<b>Medicon</b>	15.76.14
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Artery Forceps</a>
<b>KLS Martin</b>	13-248-14
<b>AUTHORS</b>	<a href="#">SPENCER</a> , <a href="#">WELLS</a>
<b>Medicon</b>	15.76.14

## SPECIFICATIONS

[serrated](#), [straight](#)