

SPENCER-WELLS HAEMOSTATIC FORCEPS

straight; serrated

15.5cm (6.2")



SKU: 040-290-15

Categories: [Hemostats](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps
Aesculap	BH334
KLS Martin	13-248-15
AUTHORS	SPENCER , WELLS
Medicon	15.76.16
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps
Aesculap	BH334
KLS Martin	13-248-15

AUTHORS

[SPENCER](#), [WELLS](#)

Medicon

15.76.16

SPECIFICATIONS

[serrated](#), [straight](#)