

SPENCER-WELLS HAEMOSTATIC FORCEPS

straight; serrated

18.0cm (7")



SKU: 040-290-18

Categories: [Hemostats](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps
Aesculap	BH336
KLS Martin	13-248-18
AUTHORS	SPENCER, WELLS
Medicon	15.76.18
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps
Aesculap	BH336
KLS Martin	13-248-18

AUTHORS

[SPENCER](#), [WELLS](#)

Medicon

15.76.18

SPECIFICATIONS

[serrated](#), [straight](#)