

## SPENCER-WELLS HAEMOSTATIC FORCEPS

straight; serrated

20.0cm (8")



**SKU:** 040-290-20

**Categories:** [Hemostats](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Artery Forceps</a>
<b>Aesculap</b>	BH338
<b>KLS Martin</b>	13-248-20
<b>AUTHORS</b>	<a href="#">SPENCER, WELLS</a>
<b>Medicon</b>	15.76.20
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Artery Forceps</a>
<b>Aesculap</b>	BH338
<b>KLS Martin</b>	13-248-20

## **AUTHORS**

[SPENCER](#), [WELLS](#)

## **Medicon**

15.76.20

## **SPECIFICATIONS**

[serrated](#), [straight](#)